

# COLLECTION STANDARD OPERATING PROCEDURE

## DISSOLVED ORGANIC CARBON (DOC) & CARBON DISSOLVED ORGANIC MATTER (CDOM)

### I. Sample collection

---

- Wear protective gloves to avoid contamination
- Samples are collected from the appropriate Niskin bottles into 250 mL amber glass bottles pre-rinsed 3 times with the sample itself

Warning: leave some space before closing the bottle;

- Samples are immediately stored at -20 °C until analysis.

### II. List of material, equipment and supplies

---

- Amber glass bottles
- Refrigerated box for transportation to the lab